

LIPOSUCTION CANNULAS



MERCEDES

Three openings in a circumferential pattern.



TOLEDO V

One hole.



STANDARD ONE HOLE

One hole near the tip.



COBRA

Two side holes near tip and one above.



INFILTRATION

10 holes in spiral pattern.

LIPOSUCTION CANNULAS



ACCELERATOR

Three openings in triangle pattern.
One center hole and two lateral holes.



CANE

Two elongated holes.



STANDARD TWO HOLE

Two hole near the tip.
In linear pattern.



SATTLER

24 holes. Eight holes on each side.



INFILTRATION

8 holes.

LIPOSUCTION CANNULAS



TWO PORT RADIAL

Two lateral holes.



CANDY CANE

Three elongated holes.



STANDARD THREE HOLES

Three holes in linear pattern.



TORNADO

Seven small holes and one hole on the tip.



COLEMAN (CURVED)

One hole near tip.



LIPOSUCTION CANNULAS

CANNULA SIZES

1.2 mm	•	MICRO CANNULA SIZES
1.4 mm	•	
1.6 mm	•	
2.0 mm	•	
2.7 mm	•	
3.0 mm	•	

3.5 mm	•	STANDARD CANNULA SIZES
4.0 mm	•	
5.0 mm	•	
6.0 mm	•	
7.0 mm	•	
8.0 mm	•	
9.0 mm	•	
10 mm	•	


CANNULA TIP STYLES

Blunt	
Bullet	
Spatula	

Note: For custom made cannulas you can send us drawings / specifications at: info@ammadsurgicals.com

CURVE PREFERENCE

- a. Convex Curve
- b. Straight
- c. Concave Curve
- d. "S" TiP



CANNULA HANDLES, HUBS & FAT TRANSFER



LUER / LUER

Female luer on top For Cannula connection & Male luer connection for syringe. Dia. 16mm, Length 5.5" Available in stainless steel & Aluminum.



LUER

Female luer on top For Cannula connection. Aluminum for light weight. Colour Option Available.



FAT TRANSFER

Male luer lock on both sides. Length: 6cm



LUER LOCK HUB (I)

Hexagonal luer lock hub for luer lock handle & syringe connection.



LUER LOCK HUB (II)

Hexagonal luer lock hub for luer lock handle & syringe connection.



LUER LOCK HUB (III)

Round, Knurled luer lock hub for luer lock handle & syringe connection. Available in aluminum.



Silver



Black



Green



Gold



Blue




Red

ALUMINUM HANDLES

Anodised
Knurled
Diameter: 25mm
Length: 7.5"

Disposable cannulas are packed in double peel away pouches sterilized. (EO) As shown.



DISPOSABLE FAT TRANSFER

Male luer lock on both sides. Pack of 10 pcs.



DISPOSABLE LUER LOCK

Hexagonal luer lock hub for luer lock handle & syringe connection.



DISPOSABLE

For single piece cannulas. Dia. 25mm, Length 8"

