

PAL Cannulas

Single-Port Cannulas

- 90° of extraction
- Non-aggressive
- Good for sensitive areas such as the neck and back
- Ports can be directed away from the skin to prevent surface irregularities



Tip Style	Multi-Use	Single-Use	Diameter	Length	Port Length
Single Port	—	PAL-300LL	3mm	30cm	8mm
	—	PAL-300LS	3mm	22cm	8mm
	—	PAL-300SS	3mm	22cm	6mm
	—	PAL-240LS	2.4mm	22cm	7.5mm
	—	PAL-240SS	2.4mm	22cm	5mm
	—	PAL-240LM	2.4mm	15cm	7.5mm

Mercedes Cannulas

- 360° of extraction
- Moderately aggressive
- The most popular cannula style



Mercedes	PAL-R504LL	PAL-504LL	5mm	30cm	10mm
	PAL-R504LS	PAL-504LS	5mm	22cm	10mm
	PAL-R404LL	PAL-404LL	4mm	30cm	10mm
	PAL-R404LS	PAL-404LS	4mm	22cm	10mm
	PAL-R304LL	PAL-304LL	3mm	30cm	8mm
	PAL-R304LS	PAL-304LS	3mm	22cm	8mm
	PAL-R304LM	PAL-304LM	3mm	15cm	8mm
	PAL-R244LM	PAL-244LM	2.4mm	15cm	7.5mm
	PAL-R404XL	—	4mm	40cm	10mm

Double Mercedes Cannulas

- 360° of extraction
- Aggressive
- Good for fibrous tissue and large-volume extraction



Double Mercedes	PAL-R505LL	PAL-505LL	5mm	30cm	10mm
	PAL-R505LS	PAL-505LS	5mm	22cm	10mm
	PAL-R405LL	PAL-405LL	4mm	30cm	10mm
	PAL-R405LS	PAL-405LS	4mm	22cm	10mm

Mirrored Tri-Port II Cannulas

- 270° of extraction
- Aggressive
- Ports can be directed away from the skin to prevent surface irregularities



Mirrored Tri-Port II	PAL-R506LL	PAL-506LL	5mm	30cm	10mm
	PAL-R506LS	PAL-506LS	5mm	22cm	10mm
	PAL-R406LL	PAL-406LL	4mm	30cm	10mm
	PAL-R406LS	PAL-406LS	4mm	22cm	10mm

Tri-Port II Cannulas

- 270° of extraction
- Moderately aggressive
- Ports can be directed away from the skin to prevent surface irregularities



Tri-Port II	PAL-R402LL	PAL-402LL	4mm	30cm	10mm
	PAL-R402LS	PAL-402LS	4mm	22cm	10mm
	PAL-R302LL	PAL-302LL	3mm	30cm	8mm
	PAL-R302LS	PAL-302LS	3mm	22cm	8mm

Tri-Port III Cannulas

- 90° of extraction
- Moderately aggressive
- Ports can be directed away from the skin to prevent surface irregularities



Tri-Port III	PAL-R403LL	PAL-403LL	4mm	30cm	10mm
	PAL-R403LS	PAL-403LS	4mm	22cm	10mm
	PAL-R303LL	PAL-303LL	3mm	30cm	8mm
	PAL-R303LS	PAL-303LS	3mm	22cm	8mm

PAL Cannulas

HARVESTING

Multi-Hole Cannulas

- 360° of extraction
- Non-aggressive
- Produces small lobules for small-volume grafting
- Also ideal for fine sculpting and feathering
- Excellent for removing fat after laser or ultrasonic procedures

Del Vecchio Track-12 Cannulas

- 270° of extraction
- Aggressive
- Designed for harvesting small lobules for grafting

SPECIALTY

Helixed Tri-Port III Cannulas

- 270° of extraction
- Aggressive
- Designed for high-fibrous areas
- Ideal for gynecomastia cases

Spatula Cannulas

- 90° of extraction
- Short, non-aggressive cannula with a blunt tip
- Ideal for delicate areas such as the neck and back

Flared Mercedes Cannulas

- 360° of extraction
- Aggressive
- Ideal for debulking, dissection, and pre-tunneling
- Good for breaking up fat without suction, prior to sculpting and extraction

Bent Cannulas

- Effective for pre-tunneling
- Good for debulking and extraction

Turbo Cannulas

- Allows air to vent through the tube during liposuction for better flow rates
- Not recommended for harvesting fat to be re-injected

HIGH DEFINITION

HD Cannulas

- 270° of extraction
- Moderately aggressive
- Ideal for sculpting of superficial fat layer

Tip Style	Multi-Use	Single-Use	Diameter	Length	Port Length
Multi-Hole	PAL-R30MLL	PAL-30MLL	3mm	30cm	1.9mm
	PAL-R30MLS	PAL-30MLS	3mm	22cm	1.9mm
	PAL-R30MLM	PAL-30MLM	3mm	15cm	1.9mm
Track-12	—	PAL-4012LL	4mm	30cm	6mm
	—	PAL-3012LS	3mm	22cm	6mm
	PAL-R4012XL	—	4mm	40cm	6mm
Helixed Tri-Port III	—	PAL-507LL	5mm	30cm	10mm
	—	PAL-507LS	5mm	22cm	10mm
	PAL-R407LL	PAL-407LL	4mm	30cm	10mm
	PAL-R407LS	PAL-407LS	4mm	22cm	10mm
Spatula (Three-Hole)	—	PAL-30S3	3mm	15cm	1.5mm
Spatula (Two-Hole)	—	PAL-30S2	3mm	15cm	1.5mm
Spatula (Single-Hole)	—	PAL-24S1	2.4mm	15cm	3mm
Flared Mercedes	PAL-R508LL	PAL-508LL	5mm	30cm	10mm
	PAL-R508LS	PAL-508LS	5mm	22cm	10mm
	PAL-R408LL	PAL-408LL	4mm	30cm	9mm
	PAL-R408LS	PAL-408LS	4mm	22cm	9mm
	PAL-R308LL	PAL-308LL	3mm	30cm	8mm
	PAL-R308LS	PAL-308LS	3mm	22cm	8mm
Bent Flared Mercedes	PAL-R508LLB	PAL-508LLB	5mm	30cm	10mm
	PAL-R408LLB	PAL-408LLB	4mm	30cm	9mm
Bent Mercedes	PAL-R404LLB	PAL-404LLB	4mm	30cm	10mm
	PAL-R304LLB	PAL-304LLB	3mm	30cm	8mm
Double Mercedes	—	PAL-505LLT	5mm	30cm	10mm
	—	PAL-405LLT	4mm	30cm	10mm
Helixed Tri-Port III	—	PAL-507LLT	5mm	30cm	10mm
Mercedes	—	PAL-404LLT	4mm	30cm	10mm
	—	PAL-404LST	4mm	22cm	10mm
Left Bent HD	—	PAL-40HLLB	4mm	30cm	10mm
	—	PAL-50HLLB	5mm	30cm	10mm
Straight HD	—	PAL-40HLL	4mm	30cm	10mm
	—	PAL-50HLL	5mm	30cm	10mm
Right Bent HD	—	PAL-40HLRB	4mm	30cm	10mm
	—	PAL-50HLRB	5mm	30cm	10mm