

PAL Cannulas

Single-Port Cannulas



- 90° of extraction
- Non-aggressive
- Good for sensitive areas such as the neck and back
- Ports can be directed away from the skin to prevent surface irregularities

Tip Style	Multi-Use	Single-Use	Diameter	Length	Port Length
Single Port	—	PAL-300LL	3mm	30cm	8mm
	—	PAL-300LS	3mm	22cm	8mm
	—	PAL-300SS	3mm	22cm	6mm
	—	PAL-240LS	2.4mm	22cm	7.5mm
	—	PAL-240SS	2.4mm	22cm	5mm
	—	PAL-240LM	2.4mm	15cm	7.5mm

Mercedes Cannulas



- 360° of extraction
- Moderately aggressive
- The most popular cannula style

Mercedes	PAL-R504LL	PAL-504LL	5mm	30cm	10mm
	PAL-R504LS	PAL-504LS	5mm	22cm	10mm
	PAL-R404LL	PAL-404LL	4mm	30cm	10mm
	PAL-R404LS	PAL-404LS	4mm	22cm	10mm
	PAL-R304LL	PAL-304LL	3mm	30cm	8mm
	PAL-R304LS	PAL-304LS	3mm	22cm	8mm
	PAL-R304LM	PAL-304LM	3mm	15cm	8mm
	PAL-R244LM	PAL-244LM	2.4mm	15cm	7.5mm
	PAL-R404XL	—	4mm	40cm	10mm

Double Mercedes Cannulas



- 360° of extraction
- Aggressive
- Good for fibrous tissue and large-volume extraction

Double Mercedes	PAL-R505LL	PAL-505LL	5mm	30cm	10mm
	PAL-R505LS	PAL-505LS	5mm	22cm	10mm
	PAL-R405LL	PAL-405LL	4mm	30cm	10mm
	PAL-R405LS	PAL-405LS	4mm	22cm	10mm

Mirrored Tri-Port II Cannulas



- 270° of extraction
- Aggressive
- Ports can be directed away from the skin to prevent surface irregularities

Mirrored Tri-Port II	PAL-R506LL	PAL-506LL	5mm	30cm	10mm
	PAL-R506LS	PAL-506LS	5mm	22cm	10mm
	PAL-R406LL	PAL-406LL	4mm	30cm	10mm
	PAL-R406LS	PAL-406LS	4mm	22cm	10mm

Tri-Port II Cannulas



- 270° of extraction
- Moderately aggressive
- Ports can be directed away from the skin to prevent surface irregularities

Tri-Port II	PAL-R402LL	PAL-402LL	4mm	30cm	10mm
	PAL-R402LS	PAL-402LS	4mm	22cm	10mm
	PAL-R302LL	PAL-302LL	3mm	30cm	8mm
	PAL-R302LS	PAL-302LS	3mm	22cm	8mm

Tri-Port III Cannulas



- 90° of extraction
- Moderately aggressive
- Ports can be directed away from the skin to prevent surface irregularities

Tri-Port III	PAL-R403LL	PAL-403LL	4mm	30cm	10mm
	PAL-R403LS	PAL-403LS	4mm	22cm	10mm
	PAL-R303LL	PAL-303LL	3mm	30cm	8mm
	PAL-R303LS	PAL-303LS	3mm	22cm	8mm

PAL Cannulas		Tip Style	Multi-Use	Single-Use	Diameter	Length	Port Length
HARVESTING	Multi-Hole Cannulas  • 360° of extraction • Non-aggressive • Produces small lobules for small-volume grafting • Also ideal for fine sculpting and feathering • Excellent for removing fat after laser or ultrasonic procedures	Multi-Hole	PAL-R30MLL	PAL-30MLL	3mm	30cm	1.9mm
			PAL-R30MLS	PAL-30MLS	3mm	22cm	1.9mm
			PAL-R30MLM	PAL-30MLM	3mm	15cm	1.9mm
	Del Vecchio Track-12 Cannulas  • 270° of extraction • Aggressive • Designed for harvesting small lobules for grafting	Track-12	—	PAL-4012LL	4mm	30cm	6mm
			—	PAL-3012LS	3mm	22cm	6mm
			PAL-R4012XL	—	4mm	40cm	6mm
SPECIALTY	Helixed Tri-Port III Cannulas  • 270° of extraction • Aggressive • Designed for high-fibrous areas • Ideal for gynecomastia cases	Helixed Tri-Port III	—	PAL-507LL	5mm	30cm	10mm
			—	PAL-507LS	5mm	22cm	10mm
			PAL-R407LL	PAL-407LL	4mm	30cm	10mm
			PAL-R407LS	PAL-407LS	4mm	22cm	10mm
	Spatula Cannulas  • 90° of extraction • Short, non-aggressive cannula with a blunt tip • Ideal for delicate areas such as the neck and back	Spatula (Three-Hole)	—	PAL-30S3	3mm	15cm	1.5mm
		Spatula (Two-Hole)	—	PAL-30S2	3mm	15cm	1.5mm
		Spatula (Single-Hole)	—	PAL-24S1	2.4mm	15cm	3mm
	Flared Mercedes Cannulas  • 360° of extraction • Aggressive • Ideal for debulking, dissection, and pre-tunneling • Good for breaking up fat without suction, prior to sculpting and extraction	Flared Mercedes	PAL-R508LL	PAL-508LL	5mm	30cm	10mm
			PAL-R508LS	PAL-508LS	5mm	22cm	10mm
			PAL-R408LL	PAL-408LL	4mm	30cm	9mm
			PAL-R408LS	PAL-408LS	4mm	22cm	9mm
			PAL-R308LL	PAL-308LL	3mm	30cm	8mm
			PAL-R308LS	PAL-308LS	3mm	22cm	8mm
	Bent Cannulas  • Effective for pre-tunneling • Good for debulking and extraction	Bent Flared Mercedes	PAL-R508LLB	PAL-508LLB	5mm	30cm	10mm
			PAL-R408LLB	PAL-408LLB	4mm	30cm	9mm
		Bent Mercedes	PAL-R404LLB	PAL-404LLB	4mm	30cm	10mm
			PAL-R304LLB	PAL-304LLB	3mm	30cm	8mm
	Turbo Cannulas  • Allows air to vent through the tube during liposuction for better flow rates • Not recommended for harvesting fat to be re-injected	Double Mercedes	—	PAL-505LLT	5mm	30cm	10mm
			—	PAL-405LLT	4mm	30cm	10mm
		Helixed Tri-Port III	—	PAL-507LLT	5mm	30cm	10mm
		Mercedes	—	PAL-404LLT	4mm	30cm	10mm
			—	PAL-404LST	4mm	22cm	10mm
	HD Cannulas  • 270° of extraction • Moderately aggressive • Ideal for sculpting of superficial fat layer	Left Bent HD	—	PAL-40HLLB	4mm	30cm	10mm
			—	PAL-50HLLB	5mm	30cm	10mm
		Straight HD	—	PAL-40HLL	4mm	30cm	10mm
			—	PAL-50HLL	5mm	30cm	10mm
		Right Bent HD	—	PAL-40HLRB	4mm	30cm	10mm
			—	PAL-50HLRB	5mm	30cm	10mm